

04

Histology



Specimen Bag



Coplin Jars



Paraffin Wax



▣ **Biohazard specimen bags** protect healthcare workers from specimen leaks while keeping its contents free from contaminants such as dirt and or moisture. Each bag is imprinted with the universal biohazard symbol.

Types

- Biohazard Specimen Bags, 13cmx25cm with pocket

▣ **Coplin jars** are wide-mouth glass jars, usually with vertically grooved interior walls, used for the storage or staining of slides containing blood smears or tissue sections.

Types

- Coplin Staining Jar, (Round), Glass, 5 Internal Slots, Holds up to 10 1"x3" slides
- Coplin Staining Jar (Square), Glass, 9 Internal Slots, Holds up to 10 1"x3" slides

▣ **Paraffin wax**, colourless or white, somewhat translucent hard wax consisting of a mixture of solid straight-chain hydrocarbons ranging in melting point from about 54°C to 56°C.

Tissue Cassette



Tissue cassettes are suitable for holding & identifying tissue samples in processing, embedding, and sectioning procedures.

Types



- Tissue Cassette Stripe Holes (White, Blue, Green, and Pink)

Cervix Examination Brush



Cervical Rambrush




✦ **Cervical Brush** provides a clear and accurate sample for the early detection of cervical cancer. Cells from the entire cervical area may be collected to produce good, representative samples.

Types

- Broom type, Cyto Brush, Pap Smear Kit/Gyne Kit

Embedding Rings



Histology Base Molds



✦ **Embedding Rings** are used in positioning histological tissues accurately in base molds.

Available colors:




✦ **Histology Base Molds** enables tissue to be positioned and embedded in wax when using routine or biopsy cassettes.

Stainless Steel

- 15 x 15 x 6 mm
- 24 x 24 x 6 mm
- 30 x 24 x 6 mm
- 37 x 24 x 6 mm

PVC Material

- 17 x 17 x 5 mm
- 25 x 25 x 5 mm
- 30 x 20 x 12 mm
- 25 x 40 x 5 mm
- 24 x 40 x 4 mm



Staining Jar and Rack

- Staining jars and staining dishes are used to stain cells and tissues on microscope slides. Without proper staining, most cells are just a clear substance.

Types: Staining Jars

1. Staining Jar, Glass, Holds a 20-Slide Rack, w/ Lid, Papanicolau Type
2. Staining Jar, Glass, Holds a 30-Slide Rack, w/ Lid
3. Staining Jar, Glass, Holds a 60-Slide Rack, w/ Lid
4. Staining Jar, Plastic, Microwave Safe, Holds a 24-Slide Rack, with lid

Types: Staining Racks

5. Staining Rack, Papanicolau Type, Stainless Steel, Holds up to 20 1"x3" slides
6. Staining Rack, Stainless Steel, Holds up to 30 1"x3" slides
7. Staining Rack, Stainless Steel, Holds up to 60 1"x3" slides
8. Staining Rack, Plastic, Microwave Safe, Holds up to 24 1" x 3" Slides