

02


Reagents & Rapid Tests



Clinical Chemistry Reagents



SGPT Kinetic 1x50ml



SGOT Kinetic 1x50ml



Cholesterol 1x200ml



Total Protein Test Kit 1x200ml



HDL Cholesterol 1x200ml



Triglyceride 1x200ml



Urea 2x100ml



Uric Acid 1x200ml



Creatinine (Kinetic) 2x100ml



Direct Bilirubin Test Kit 1x200ml



Total Bilirubin Test Kit 1x200ml



Glucose 1x200ml



Biological Stains



Gram Stain 4x250ml



Rapi Stain 3x500ml



Giemsa Stain 1x1000ml



Serology Blood Grouping Reagents



Anti B Sera (1x10ml)



Anti A Sera (1x10ml)



Anti D(Rh) Sera (1x10ml)

RAPID TESTS

Designed for use where a preliminary screening test result is required.




Laboquick Pregnancy Test is an immunochromatographic based one step quick test that is aimed at qualitative detection of human chorionic gonadotropin hormone (hCG) in urine/serum.



Laboquick Fecal Occult Blood Test is a lab test used to check stool samples for hidden (occult) blood. The occult blood in the stool may indicate colon cancer or polyps in the colon or rectum - though not all cancer or polyps bleed.



Laboquick Troponin 1 measures the levels of Troponin 1 proteins in the blood. These proteins are released when the heart muscle has been damaged, such as occurs with a heart attack.



Wondfo Drug of Abuse Test Cards are intended for the qualitative detection of DOA markers and their principal metabolites in human urine at specified cut-off level.



Standard Q Rapid Tests

Made in Korea



DENGUE DUO

STANDARD Q Dengue Duo test is an immunochromatographic assay for the detection of NS1 antigen and Dengue virus-specific IgM and IgG antibodies in human serum, plasma and whole blood.

Specifications:

Specimen: Serum, Plasma, Whole Blood

Sensitivity: NS1 (92.9%, 184/198), IgM (97.5%, 77/79), IgG (97.2%, 140/144)

Specificity: NS1 (98.7%, 222/25), IgM (96.6%, 346/358), IgG (96.2%, 282/293)

Specimen Volume: 10µl (Serum, Plasma - Dengue IgM, IgG), 100µl (Whole blood - Dengue NS1 Ag)

Test time: 15-20mins (Inaccurate results after 20 minutes)

Storage temperature: 2-40° C/ 34-104°F



DENGUE IgM/IgG

STANDARD Q Dengue IgM/IgG Test is an immunochromatographic assay for the detection of IgM/IgG antibodies to dengue in human serum, plasma, or whole blood samples.

Specifications:

Specimen: Whole blood, Serum, Plasma

Sensitivity: IgM (97.5%, 77/79), IgG (97.2%, 140/144)

Specificity: IgM (96.6%, 346/358), IgG (96.2%, 282/293)

Test time: 15-20mins (Inaccurate results after 20 minutes)

Storage temperature: 2-40° C/ 34-104°F

Specimen Volume: 10µl



DENGUE NS1 Ag

STANDARD Q Dengue NS1 Ag Test is an immunochromatographic assay for the detection of NS1 antigens in human serum, plasma, or whole blood samples.

Specifications:

Specimen: Whole blood, Serum, Plasma

Sensitivity: 92.9% (184/198)

Specificity: 98.7% (222/225)

Test time: 15-20mins (Inaccurate results after 20 minutes)

Storage temperature: 2-40° C/ 34-104°F



Covid-19 Ag (Nasopharyngeal)

STANDARD Q COVID-19 Ag Test is a rapid chromatographic immunoassay for the qualitative detection of nucleocapsid antigens to SARS-CoV-2 present in human nasopharynx.

Specifications:

Specimen: Nasopharyngeal

Sensitivity: Ct ≤ 25 : 97.14%(68/70), Ct ≤ 33 : 90.71%(127/140, overall : 84.97%(130/153)

Specificity: 98.94% (1490/1506)

Test time: In 15~30 minutes

Storage temperature: 2-30° C/ 36-86°F

Specimen Volume: 3 drop



Covid-19 Ag (Nasal)

STANDARD Q COVID-19 Ag Test is a rapid chromatographic immunoassay for the qualitative detection of nucleocapsid antigens to SARS-CoV-2 present in human nasal or nasopharyngeal specimens.

Specifications:

Specimen: Nasal & Nasopharyngeal

Sensitivity: 97.12% (101/104)

Specificity: 100% (399/399)

Test time: In 15~30 minutes

Storage temperature: 2-30° C/ 36-86°F

Specimen Volume: 3 drop

Glucose Drink



The **Oral Glucose Tolerance Test (OGTT)** was the gold standard for making the diagnosis of type 2 diabetes. It is still commonly used during pregnancy for diagnosing gestational diabetes. With an oral glucose tolerance test, the person fasts overnight (at least 8 hours, but not more than 16 hours).

Product Features:

- BVO free
- Certified Kosher
- Gluten and dairy free
- Shatter-proof plastic bottle
- Easy-to-read volume indicator
- Sizes: 50g, 75g, 100g