

10

General Laboratory Products



Filter Paper

✦ Used for routine applications with medium retention and flow rate.

Size

600mm x
600mm

Weight

80
grams/sqm

Format


Ashless

Digital Timer


✦ This timer may be used for various timing purposes in the laboratory.

Specifications

- Battery: AAA alkaline battery
- Single Channel, 60 sec Sound Alarm
- Can be set from 1 sec to 99 min 59 sec
- Pocket Clip, Stand Allow Multiple Use
- Dimension: 79 x 76 x 18mm



Refrigerator Thermometer



Alcohol Filled Thermometer



Laboratory Thermometer

✦ An instrument for determining temperatures.

Alcohol Filled

- With Red Spirit- Filled Celsius Scale
- Range -10°C to 110°C, 1°C Divisions

Refrigerator Thermometer

Range: 30° to 30°C/ -20° to 80° F
Diameter: Φ 52mm
Material: Stainless Steel



Lens Cleaning Tissue

✦ Can be used to clean lenses and other glass objects without scratching the surface.

Brand: Sterilab

Weight: 11gr/sq.meter

Dimension: 150 mm

Packaging: 100 sheets/pad



Parafilm

✦ Flexible, thermoplastic used to prevent adhering and for a wide variety of clinical applications

Brand: Parafilm M

Size: 4 inches × 125 feet

Temperature Limit: -45 to +50 °C



pH Indicator Paper

✦ It is used to determine if a solution is acidic, basic or neutral. This is determined by dipping part of the paper into a solution of interest and watching the color change.

Brand: Sterilab

Size: 82 gram/sqm

Range: 1-14



Needle Destroyer

✦ Helps to destroy needles and syringes immediately after use.

Types:

- With cutter, Round
- With cutter, Square



Sharps Container

✦ Sharps containers are commonly used across the health care industry to dispose of needles, syringes, and other materials that can pierce or cut the skin.

Sizes:

- 5 liters
- 7 liters



Test Tube Brushes

✦ Test tube brushes are used for cleaning test tubes. These brushes are very durable and have hog bristles and a wire body.

Types:

- Small
- Medium
- Large

Partner Brands

